

# Autoclave/Steam Indicator Label

## PMAL7025

# DATA SHEET

### Description

Steam sterilization process indicators are Class 1 process indicators and are manufactured in compliance with ISO 11140-1:2005. The indicator will provide a distinct colour change (Red to black) after being subjected to dry steam for 10 minutes at 121°C, or 3 minutes at 134°C



### Feature

- Consistent colour change on a uniform label from red to black
- Each batch is validated in-house
- The chemistry and product ranges have been completely developed in house

### Details

12 mm Self Adhesive labels on a reel  
 \*\*The colours are for guidance only and are not authentic reproductions

### Example



### Exposed



### Shelf Life & Storage Conditions

Sterisure steam Indicator labels have a shelf life of 18 months from the date of manufacture prior to use, provided they are stored in ambient conditions below 20°C and 50% R.H. and away from direct sunlight.

### Adhesive Application Guide

Polycarbonate •	ABS •	Cleanroom Paper
Polyester •	Rubber	Carton Box •
Polyethylene •	Glass •	
Polypropylene •	Metal •	

**Format:** Self-adhesive label 12 mm diameter, 5000 per reel, 1 dot per width